

Parts List

Headband x1 Foam Pad x1 Visor Band x1

Visor Sheet x4 Tape and Adjuster x1

Assembly Instructions

<p>1. Remove adhesive backing from foam pad and attach to headband inside the guidelines marked on the headband.</p>	<p>2. Thread elastic tape through the outer slot in headband as shown. Note: Orientation of the Adjuster.</p>
<p>3. Thread elastic tape through the inner slot in the Headband as shown.</p>	<p>4. Bend the Headband (it can help to hold it between your knees) and thread the tape through the Inner and outer slots at other end of the Headband as shown</p>
<p>5. Thread the tape through the Adjuster as shown</p>	<p>6. Attach Left hand end of visor band to Left side of head band as shown. Note Letters 'L' must be aligned upright on both parts.</p>
<p>5. Rotate the visor band ~30degrees anti-clockwise past the detent to lock in place</p>	<p>6. Bend the Visor band to attach the right-hand end to the pivot on the Head band. Note the letters R must be aligned upright on both parts.</p>

<p>7. Rotate the right-hand end of the visor band ~30degree clock-wise past the detent to lock in place.</p>	<p>8. Attach the visor sheet at one side of the visor band. Ensure the visor sheet is hooked under the guide tab as shown</p>
<p>9. Wrap visor sheet around, aligning holes with guide pegs and hook under second guide tab</p>	<p>10. Use the elastic tap adjuster to achieve a comfortable fit</p>
<p>Cleaning Immerse in warm water with a mild detergent and scrub with a non-abrasive cloth or soft brush.</p>	

CAREFULLY READ THESE USER INSTRUCTIONS BEFORE USING THIS VISOR

Applicable to products: Nifco Visor 1

Manufacturer: Nifco UK Limited, Eaglescliffe, Stockton on Tees, TS16 0PS

Notified body: Module B and C2 - SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, D15YN2P. Republic of Ireland (Notified Body 2777). Declaration of conformity (www.nifcoeu.com)

These products are classed as Category III Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425 and have been shown to comply with this Regulation through limited testing against the Harmonised European Standard EN 166:2001 (see below).

This product is designed to provide some protection against contact with airborne droplets that may constitute a biological hazard. However, always remember that no item of PPE can provide full protection and care must always be taken while carrying out the risk-related activity. Never use this visor for protection against other hazards for which the product is not design for or tested against

PERFORMANCE AND LIMITATIONS OF USE – These products have been tested in accordance with EN 166: 2001 Clauses: 6 † Design; 7.1.1 † Field of vision; 7.1.2.1 † Refractive powers; 7.1.2.2 † Transmittance; 7.1.2.3 † Diffusion of light; 7.1.3 † Quality of materials; 7.1.4.2 † Increased robustness; 7.1.5.1 † Stability at elevated temperature; 7.1.7 † Resistance to ignition; 7.2.4 † Droplets and splashes of liquids.

FITTING AND SIZING – To put on and take off the product, simply place the head band over your head with the visor in front of your face and adjust the fit by tightening or loosening the elastic tape through the Adjuster. Only wear products of a suitable size that are correctly adjusted. Products which are too loose may not stay in place and if too tight will be uncomfortable. The length and tension on the elastic tape can be adjusted to fit most adults.

COMPATIBILITY – To optimise protection, it will be necessary to use these products with suitable gloves/ gowns/ masks. In this case, before carrying out the risk-related activity, consult your supplier or supervisor to ensure that all your protective products are compatible and suitable for your application. Warning - Parts of the visor that come into contact with your skin may lead to allergic reactions. In this case please discontinue use and seek medical advice.

STORAGE AND TRANSPORT – When not in use, store the product in a dark, well-ventilated area away from extremes of temperature and sunlight. Never place heavy items on top of it. If possible, avoid excessive folding and preferably store it hanging vertically. If the product is wet, allow it to dry fully before placing it into storage.

OBSOLESCENCE - This visor has a shelf life of 3 years. Service life will depend on usage conditions, but the visor should be disposed of and replaced if badly scratched or damaged in some way that prevents safe use. The visor sheet is replaceable. Replacements can be purchased at www.nifcoeu.com.

Expires: date

REPAIR – If the visor becomes damaged, it will NOT provide the optimum level of protection, and therefore should be immediately replaced. All replacements can be purchased at www.nifcoeu.com If in doubt always consult the manufacturer before attempting any type of repair.

CLEANING – Only use the following recommended cleaning treatment - Rinse with warm water and a mild detergent. NEVER use solvent based cleaning agents. New, replacement visors can be inserted into the headband if they become damaged.